



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR146-2



Product Description:

high-speed laboratory centrifuge is a routine laboratory centrifuge, widely used in research, education and production departments of biology, chemistry, medicine, etc. It is used in laboratories, scientific research, schools and other departments, and is suitable for rapid separation and synthesis of trace samples..

Large capacity refrigerated centrifuge technical performance:

- ✓ The max speed is 18500rpm, and the largest capacity of swing rotor is 4X500ml.
- ✓ It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Automatic rotor recognition make safety for multiple rotors choice. (Option)
- ✓ Up to 9 user programs storage meet more demands.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#Let'sgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR146-2



Technical Specifications :

Model	YR146-2	YR131-2
Max speed	18500r/min	5000r/min
Max RCF	23797xg	5030xg
Max capacity	4x500ml	
Speed accuracy	± 10r/min	
Time setting	1s~99min59s	
Noise	<65db(A)	
Power supply	AC110V/AC220V, 50Hz-60Hz	
Consumption	750w	
Dimension(L×W×H)	465mm×560mm×380mm	
Packing size(L×W×H)	560mm×650mm×470mm	
Net weight (No rotor)	50kg	



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges

Tabletop High Speed Centrifuge

Model YR146-2



Rotor for YR146-2

 <p>No.1</p>	<p>Capacity: 12x1.5ml Max Speed:18500r/min Max RCF: 23797xg Adaptor: 0.2ml、0.5ml</p>	 <p>No.2</p>	<p>Capacity: 10x5ml Max Speed: 13000r/min Max RCF: 11400xg Adaptor: 0.5ml、1.5ml</p>
 <p>No.2</p>	<p>Capacity: 12x5ml Max Speed: 13000r/min Max RCF: 11400xg</p>	 <p>No.3</p>	<p>Capacity: 12x10ml Max Speed:12000r/min Max RCF: 14800xg Adaptor: 1.5ml、5ml</p>
 <p>No.3</p>	<p>Capacity: 8x15ml Max Speed:12000r/min Max RCF: 14800xg Adaptor: 1.5ml、5ml</p>	 <p>No.4</p>	<p>Capacity:6x50ml Max Speed: 13000r/min Max RCF: 17940xg</p>
 <p>No.4</p>	<p>Capacity: 4x100ml Max Speed: 12000r/min Max RCF: 13000xg Adaptor: 10ml、20ml</p>	 <p>No.7</p>	<p>Capacity: 24x1.5ml/2.2ml Max Speed: 13200r/min Max RCF: 16110xg Adaptor: 0.2ml、0.5ml</p>
 <p>No.8</p>	<p>Capacity:36x0.5ml Max Speed: 16000r/min Max RCF: 23669xg</p>	 <p>No.8</p>	<p>Capacity:48x0.5ml Max Speed:16000r/min Max RCF: 23669xg</p>
 <p>No.9</p>	<p>Capacity:48x1.5ml Max Speed: 13000r/min Max RCF: 18505xg</p>		



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.






All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges Tabletop High Speed Centrifuge Model YR146-2



Rotor for YR146-2 and YR131-2

 <p>No.1</p>	<p>Capacity:4 x250ml Max Speed:4000r/min Max RCF:3200 xg</p>	<p>250 ml Bucket 250ml PE Bottle</p>	<p>Adaptor:4x9x1.5ml、4 x7x5ml、4 x9 x5ml/7ml (Blood Vacuum Tube)、4 x 7 x 10ml (Internal Thread)、4 x 5x 15ml、4 x4 x20ml、4 x50ml、4 x100ml、4 x250ml</p>
 <p>No.2</p>	<p>Capacity:4 x500ml Max Speed:4000r/min Max RCF:3040 xg</p>	<p>4 x14 x5ml/7ml Aerosol Cover (Option) 4 x24 x5ml Blood Vacuum Tube 4 x19 x5ml/7ml Aerosol Cover (Option) 4 x37 x7ml G RIA Tube</p>	<p>Square Bucket Adaptor:4 x25 x1.5ml、4 x28 x5ml (Blood Vacuum Tube)、4 x 20 x 10ml (Internal Thread)、4 x 10x 15ml、4 x4 x50ml、4 x250ml 500 ml Bucket Adaptor: 4 x19 x5ml/7ml、4 x9x15ml、4 x 7 x 20ml、4 x3x50ml、4 x100ml、4 x250ml</p>
 <p>No.2</p>	<p>Capacity: 8 x50ml Max Speed:5000r/min Max RCF:5030xg</p>	<p>4 x 2 x 50ml</p>	<p>Adaptor: 8 x 2 x 1.5ml、8 x 5ml、8 x 5ml/7ml、8 x 10ml、8 x 15ml、8 x 20ml</p>



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.







All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S




KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR146-2



 <p>No.3</p>	<p>Max Speed: 4000r/min Max RCF:3020 xg</p>	<p>Adaptor: 4x2x1.5ml、4x4x1.5ml、4x6x1.5ml、 8x4x1.5ml、4x2x5ml/7ml、4x4x5ml/7ml、4 x6x5ml/7ml、8x4x5ml/7ml、4 x 2 x 10ml、4 x 4 x 10ml、4 x 6 x 10ml、8x4x10ml、4 x 2 x15ml、4 x2 x20ml、4 x4 x15ml、4 x6 x15ml、8x4x15ml、4 x2 x50ml</p>	
 <p>No.3</p>	<p>Max Speed: 4000r/min Max RCF:3200 xg</p>	<p>Adaptor:4 x12 x7ml、4 x6 x15ml、4 x9 x15ml、4 x50ml、4 x2 x50ml、4 x100ml</p>	
 <p>No.3</p>	<p>Capacity:4 x 18 x 7ml Max Speed:4000r/min Max RCF:3200xg</p>	 <p>4x18x7ml</p>	



#Letsgivemore 
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.




All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S




KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR146-2



 <p>No.5</p>	<p>Capacity: 4 x250ml Max Speed: 4000r/min Max RCF:3310xg</p>	<p>Adaptor: 4 x4x15ml、 4 x50ml、 4 x250ml</p>  <p>500ml</p>  <p>4x500ml</p>  <p>4x250ml</p>  <p>4x4x50ml</p>  <p>4x7x15ml</p>	
 <p>No.5</p>	<p>Capacity:4 x500ml Max Speed:4000r/min Max RCF:3310xg</p>	<p>500ml</p>	<p>500 ml Bucket Adaptor: 4 x7x15ml、 4 x4x50ml 、 4 x250ml、 4 x500ml</p>  <p>500ml</p>  <p>4x50ml</p>  <p>4x4x15ml</p>



#Letsgivemore 
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>

<https://www.kalstein.fr/>

