



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.




Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Tabletop High-Speed Refrigerated Centrifuge****Model YR0125****Product Description:**

The main working principle of high speed freezing centrifuge: high speed freezing centrifuge is the powerful centrifugal force generated by high-speed rotation of the high speed freezing centrifuge frozen centrifuge rotor to accelerate the sedimentation speed of particles in the liquid, and separate the sedimentation coefficient and buoyant density of the sample.

**Technical performance:**

- ✓ It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Imported compressor unit with Non-CFC refrigerant R134a.
- ✓ It has self-locking, over-speed detection, overtemperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.


**#Let'sgive more** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges**Tabletop High-Speed Refrigerated Centrifuge****Model YR0125****Technical Specifications :**

Model	YR0125	YR127	YR128
Max speed	20000r/min	18000r/min	16000r/min
Max RCF	27800xg	23755xg	17800xg
Max capacity	4x100ml	4x100ml	6x50ml
Speed accuracy	± 30r/min		
Time setting	1min~99h 59min		
Temperature setting	-20 °C ~ +40 °C		
Temperature accuracy	±1 °C		
Compressor unit	French Tecumseh Compressor Unit R134a		
Noise	<65db(A)		
Power supply	AC110V/AC220V, 50Hz-60Hz		
Consumption	1.2kw		
Dimension(L×W×H)	600mm×470mm×330mm		
Packing size(L×W×H)	720mm×640mm×430mm		
Net weight (No rotor)	50kg		

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Tabletop High-Speed Refrigerated Centrifuge
Model YR0125



CENTRIFUGES

 	<p>Capacity: 12x1.5ml/2.2ml Max Speed: 20000/18000/16000(r/min) Max RCF: 27800xg/22530xg/17800xg Adaptor: 0.2ml、 0.5ml</p>	 	<p>Capacity: 24x1.5ml/2.2ml Max Speed: 16000/16000/13200(r/min) Max RCF: 23755xg/23755xg/16110xg Adaptor: 0.2ml、 0.5ml</p>
 	<p>Capacity: 10x5ml Max Speed: 16000/13000/13000(r/min) Max RCF: 17280xg/11400xg/11400xg Adaptor: 0.5ml、 1.5ml</p>	 	<p>Capacity: 12x5ml Max Speed: 16000/13000/13000(r/min) Max RCF: 17280xg/11400xg/11400xg Adaptor: 0.5ml、 1.5ml</p>
 	<p>Capacity: 2 x48 Max Speed: 3000r/min Max RCF:1400xg</p>	 	<p>Capacity: 36x0.5ml Max Speed: 13000/13000/12000(r/min) Max RCF: 16760xg/16760xg/14300xg</p>
 	<p>Capacity: 48x0.5ml Max Speed: 13000r/min 13000r/min 12000r/min Max RCF: 16760xg 16760xg 14300xg</p>	 	<p>Capacity: 12x10ml Max Speed: 13000r/min 13000r/min 12000r/min Max RCF: 16760xg 16760xg 14300xg Adaptor: 1.5ml、 5ml</p>
 	<p>Capacity: 8x15ml Max Speed: 13000r/min 13000r/min 12000r/min Max RCF: 16760xg 16760xg 14300xg Adaptor: 1.5ml、 5ml</p>	 	<p>Capacity: 6x50ml Max Speed: 12000r/min 12000r/min 11000r/min Max RCF: 13000xg 13000xg 10950xg Adaptor: 15ml、 20ml</p>
 	<p>Capacity: 4x100ml Max Speed: 12000r/min Max RCF: 13000xg Adaptor: 10ml、 20ml Only for TGL- 20M&TGL-18M</p>	 	<p>Capacity: 4x10ml Max Speed: 6000r/min Max RCF: 3200xg</p>
 	<p>Capacity: 48x1.5ml Max Speed: 15000/14000(r/min) Max RCF: 24637xg/21461xg Adaptor: 0.2ml、 0.5ml</p>		



#Let'sgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
https://kalstein.eu
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

