



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Tabletop High-Speed Refrigerated Centrifuge****Model YR126-1****Product Description:**

The high-speed benchtop centrifuge mainly extracts the required components from the liquid mixture. According to the density of each substance, after high-speed rotation, the dense liquid sinks to the bottom layer, and the small density liquid floats on it, thus stratifying the liquid. Refine the pure things I need. It is used in many fields such as medicine and chemical industry.

**Technical performance:**

- ✓ The max speed is 21000rpm, and the largest capacity of swing rotor is 4X750ml.
- ✓ It has micro-computer control, electronic lock and AC variable conversion motor which delivers high torque.
- ✓ LCD touch screen make visualized operation conveniently.
- ✓ Imported compressor unit with Non-CFC refrigerant R404a.
- ✓ Automatic rotor recognition make safety for multiple rotors choice. (Option)
- ✓ Up to 6 user programs storage in the surface is selected directly.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges**Tabletop High-Speed Refrigerated Centrifuge****Model YR126-1****Technical Specifications :**

Model	YR126-1	YR137-1
Max speed	21000r/min	5800r/min
Max RCF	30642xg	5060xg
Max capacity	4x750ml	
Speed accuracy	± 10r/min	
Time setting	1s~99h 59min59s	
Temperature setting	-20 °C ~ +40 °C	
Temperature accuracy	±1 °C	
Compressor unit	French Tecumseh Compressor Unit R404a	
Noise	<65db(A)	
Power supply	AC110V/AC220V, 50Hz-60Hz	
Consumption	1.5kw	
Dimension(L×W×H)	734mm×620mm×445mm	
Packing size(L×W×H)	910mm×720mm×620mm	
Net weight (No rotor)	100kg	

**#Letsgivemore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges**Tabletop High-Speed Refrigerated Centrifuge****Model YR126-1****Rotor for YR0126-1:**

 <p>No.1</p>	<p>Capacity: 12x1.5ml/2.2ml Max Speed: 21000r/min Max RCF:30642xg Adaptor: 0.2ml、0.5ml</p>	 <p>No.2</p>	<p>Capacity: 12x5ml Max Speed: 16500r/min Max RCF: 19175xg Adaptor: 0.2ml、0.5ml</p>
 <p>No.3</p>	<p>Capacity: 6x30ml Max Speed: 15000r/min Max RCF: 23120xg</p>	 <p>No.3</p>	<p>Capacity: 12x10ml Max Speed: 15000r/min Max RCF: 23120xg Adaptor: 1.5ml、5ml</p>
 <p>No.4</p>	<p>Capacity: 6x50ml Max Speed: 13000r/min Max RCF: 17940xg</p>	 <p>No.7</p>	<p>Capacity: 24x1.5ml/2.2ml Max Speed: 16500r/min Max RCF: 26050xg Adaptor: 0.2ml、0.5ml</p>
 <p>No.8</p>	<p>Capacity: 6x85ml Max Speed: 12000r/min Max RCF: 15805xg Adaptor: 5ml、15ml、20ml</p>	 <p>No.9</p>	<p>Capacity: 6x100ml Max Speed: 10000r/min Max RCF: 11960xg Adaptor: 10ml、20ml</p>
 <p>No.10</p>	<p>Capacity: 48x15ml Max Speed: 15500r/min Max RCF: 26307xg Adaptor: 0.2ml、0.5ml</p>		

**#Letsgivemore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Tabletop High-Speed Refrigerated Centrifuge

Model YR126-1



Rotor for YR0126-1 & YR137-1:

 <p>No.1</p>	<p>Capacity: 6x50ml Max Speed: 5800r/min Max RCF: 5060 xg Adaptor: 15ml, 20ml</p>	 <p>No.2</p>	<p>Capacity: 4 x2 x96 Max Speed: 4000r/min Max RCF:3040 xg Microplate Bucket 4 x2 x96 Thickness≤15.5ml x2</p>
 <p>No.3</p>	<p>Capacity: 2 x2 x96 Max Speed: 4000r/min Max RCF:2200 xg 2 x2 x96 Microplate Bucket</p>	 <p>No.7</p>	<p>Capacity:4x500ml (Conical) Max Speed:4000r/min Max RCF:3310 xg 500ml Adaptor:4 x500ml、 4 x250ml 、 4 x4 x50ml、 4 x7 x15ml</p>
 <p>No.7</p>	<p>Capacity:4 x250ml (Conical) Max Speed:4000r/min Max RCF:3310 xg 250ml Adaptor:4 x50ml、 4 x4 x15ml</p>	 <p>250ml</p>	 <p>4x50ml</p>  <p>4x4x15ml</p>



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A.S









KALSTEIN® Centrifuges

Tabletop High-Speed Refrigerated Centrifuge

Model YR126-1



 <p data-bbox="263 1022 354 1074">No.2</p>	<p>Capacity:4 x500ml Max Speed:4000r/m in Max RCF:3040 xg</p>	<p>4 x14 x5ml/7ml Aerosol Cover (Option) 4 x24 x5ml Blood Vacuum Tube 4 x19 x5ml/7ml Aerosol Cover (Option) 4 x37 x7ml G RIA Tube</p>	<p>Square Bucket Adaptor:4 x25 x1.5ml、4 x28 x5ml (Blood Vacuum Tube)、4 x 20 x 10ml (Internal Thread)、4 x 10x 15ml、 4 x4 x50ml、4 x250ml 500 ml Bucket Adaptor: 4 x19 x5ml/7ml、4 x9x15ml、4 x 7 x 20ml 、4 x3x50ml、4 x100ml 、4 x250ml</p>
	 <p data-bbox="815 1249 1252 1358">4x14x5ml/7ml Aerosol Cover (Option)</p>  <p data-bbox="859 1707 1168 1764">4x19x5ml/7ml</p>	 <p data-bbox="1306 1249 1688 1358">4x24x5ml Blood Vacuum Tube</p>  <p data-bbox="1379 1707 1597 1764">4x37x7ml</p>	



#Letsgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Tabletop High-Speed Refrigerated Centrifuge

Model YR126-1



 <p>No.4</p>	<p>Capacity:4 x750ml Max Speed:4000r/min Max RCF:3780xg</p>	<p>750ml Bucket 750ml PE Bottle</p>	<p>4 x24 x5ml/7ml (Blood Vacuum Tube)、 4x37x5ml/7ml(G RIA Tube) 4 x14 x15ml、 4 x12 x20ml、 4 x5 x50ml、 4 x3 x100ml、 4 x250ml、 4 x500ml</p>
<div style="display: flex; flex-wrap: wrap; justify-content: space-around;"> <div style="text-align: center; margin: 10px;">  <p>750 ml Bucket</p> </div> <div style="text-align: center; margin: 10px;">  <p>750ml PE Bottle</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x24x5ml/7ml Blood Vacuum Tube</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x37x5ml/7ml G RIA Tube</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x14x15ml</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x12x20ml</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x5x50ml</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x3x100ml</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x250ml</p> </div> <div style="text-align: center; margin: 10px;">  <p>4x500ml</p> </div> </div>			



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S











KALSTEIN® Centrifuges

Tabletop High-Speed Refrigerated Centrifuge

Model YR126-1



 <p>No.5</p>	<p>Capacity:6 x250ml Max Speed:4200r/min Max RCF:3500xg</p>	<p>250 ml Bucket</p>	 <p>250 ml Bucket</p>  <p>250ml PE Bottle</p>
 <p>No.6</p>  <p>No.6</p>	<p>Capacity:4 x250ml Max Speed:4300r/min Max RCF:3582xg</p>	<p>250ml PE Bottle</p>	<p>6x9x1.5ml、 6x7 x5ml、 6x9x5/7ml(Blood Vacuum Tube)、 6x7x10ml(Internal Thread)、 6x5x15ml、 6x4 x20ml、 6 x50ml、 6 x100ml</p>  <p>6x9x1.5ml</p>  <p>6x7x5ml</p>  <p>6x9x5/7ml Blood Vacuum Tube</p>  <p>6x7x10ml Internal Thread</p>  <p>6x5x15ml</p>  <p>6x4x20ml</p>  <p>6x50ml</p>  <p>6x100ml</p>



#Letsgivemore ❤️

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

