



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Tabletop Low Speed Benchtop Centrifuge****Model YR0140****Product Description:**

Low Speed Benchtop Centrifuge rpm is lower than 1000r/min, centrifugal total capacity of up to 3000ml. In clinical medicine, biochemistry, immunology, drug purification involving a wide range of applications.

Product Details:

Low Speed Benchtop Centrifuge rpm is lower than 1000r/min centrifugal total capacity of up to 3000ml. In clinical medicine biochemistry, immunology, drug purification involving a wide range o applications.

Low Speed Benchtop Centrifuge definition: Centrifuges with a speed lower than a certain level are called Low Speed Benchtop Centrifuge. The structure design is reasonable, the sanitary corners are effectively eliminated, the structural parts are smooth and smooth, the surface is polished, and the external structural parts, fasteners and drums are all rusted. Made of steel. The outer casing is of the clamshell type, which can thoroughly clean the centrifuge casing without rotating the drum sandwich space.

Technical Performance:

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy.
- ✓ Touch panel, programmable operation, operation data can be set according to demand and stored automatically.
- ✓ LED or LCD display. user-friendly interface, simple and convenient operations.
- ✓ The real-time conversation Of Speed and RCF is convenient for operation.
- ✓ It has self-locking device, lid safety device, over-speed safety device, and automatic alarm device, etc.
- ✓ Silastic sealed ring comfort to GMP.US FDA certificated.

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges**Tabletop Low Speed Benchtop Centrifuge****Model YR0140****Technical Specifications:**

Model	YR0140
Max speed	5000r/min
Max RCF	4100×g
Max capacity	15ml×16
Speed accuracy	± 30r/min
Time setting	1min~99min
Noise	<65dB(A)
Power supply	110V/60Hz 220V/50Hz
Dimension(LxWxH)	410mm×490mm×310mm
Packing size(LxWxH)	500mm×580mm×400mm
Net weight	22kg

**#Letsgivemore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S













KALSTEIN® Centrifuges

Tabletop Low Speed Benchtop Centrifuge

Model YR0140



<p>No.1</p> 	<p>Rotor Capacity: 15mlx12 Max Speed: 6000r/min Max RCF: 5120xg</p>	<p>No.1</p> 	<p>Rotor Capacity: 10mlx24 Max Speed: 5000r/min Max RCF: 4100xg</p>				
<p>No.2</p>  <p>Max Speed: 4200r/min Max RCF: 2760xg</p>	 5mlx2x12 Bucket  10mlx2x12 Bucket  15mlx12 Bucket	<p>5ml/7ml Vacuum Tube</p>			<p>Adapter: 0.5ml</p>	<p>Adapter: 1.5ml</p>	<p>Adapter: 5ml/7ml Vacuum Tube</p>
<p>No.2</p>  <p>Rotor Capacity: 15mlx16 Max Speed: 4200r/min Max RCF: 2760xg</p>	 15mlx2x8 Bucket	<p>Adapter: 1.5ml</p>		<p>Adapter: 5ml</p>			
<p>No.2</p>  <p>Rotor Capacity: 50mlx4 Max Speed: 4200r/min Max RCF: 2760xg</p>	 50mlx4 Bucket	<p>Adapter: 1.5mlx2x4</p>	<p>Adapter: 5mlx4</p>	<p>Adapter: 10mlx4</p>			
		<p>Adapter: 15mlx4</p>	<p>Adapter: 20mlx4</p>	<p>Adapter: 5ml/7mlx4 Vacuum Tube</p>			



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

