



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges

Tabletop Low Speed Centrifuge

Model YR0139



Product Description:

YR0139 Tabletop Low Speed Centrifuge is newly developed by our company on the basis of rich experience in centrifuge line. It can meet the requirement of clinical medicine..

Product Details:

YR0139 Tabletop Low Speed Centrifuge is newly developed by our company on the basis of rich experience in centrifuge line. It can meet the requirement of clinical medicine.

Low speed centrifuge definition: Centrifuges with a speed lower than a certain level are called low-speed centrifuges. The structure design is reasonable, the sanitary corners are effectively eliminated, the structural parts are smooth and smooth, the surface is polished, and the external structural parts, fasteners and drums are all rusted. Made of steel. The outer casing is of the clamshell type, which can thoroughly clean the centrifuge casing without rotating the drum sandwich space.

Technical Performance:

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy.
- ✓ Touch panel, programmable operation, operation data can be set according to demands and stored automatically.
- ✓ LED or LCD display, user-friendly interface, simple and convenient operations.
- ✓ The real-time conversation Of Speed and RCF is convenient for operation.
- ✓ It has self-locking device lid safety device, over-speed safety device and automatic alarm device, etc.
- ✓ Silastic sealed ring comfort to GMP.US FDA certificated.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges
Tabletop Low Speed Centrifuge
Model YR0139



Technical Specifications:

Model	YR0139
Max speed	4200r/min
Max RCF	2760×g
Max capacity	20ml×12
Speed accuracy	± 30r/min
Time setting	1min~90min
Noise	<65dB(A)
Power supply	110V/60Hz 220V/50Hz
Dimension(LxWxH)	410mm×490mm×310mm
Packing size(LxWxH)	500mm×580mm×400mm
Net weight	22kg



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.










All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges Tabletop Low Speed Centrifuge Model YR0139



<p>No.1</p>  <p>Optional: Aluminum Rotor Stainless Steel Tubes</p>	<p>Rotor Capacity: 10mlx18 Max Speed: 4000r/min Max RCF: 2200xg</p>	<p>Adapter: 1.5ml</p>	<p>Adapter: 5ml</p>	
<p>Adapter: 5ml/7ml Vacuum Tube</p>				
<p>No.1</p>  <p>Optional: Aluminum Rotor Stainless Steel Tubes</p>	<p>Rotor Capacity: 20mlx12 Max Speed: 4000r/min Max RCF: 2200xg</p>	<p>Adapter: 5ml</p>	<p>Adapter: 10ml</p>	
<p>Adapter: 5ml/7ml Vacuum Tube</p>				
<p>No.2</p>  <p>Max Speed: 4200r/min Max RCF: 2760xg</p>	<p> 5mlx2x8 Bucket</p> <p> 10mlx8 Bucket</p>	<p>Apply to 5ml/7ml Vacuum Tube</p>		
<p>No.2</p>  <p>Rotor Capacity: 50mlx4 Max Speed: 4200r/min Max RCF: 2760xg</p>	<p> 50mlx4 Bucket</p>	<p>Adapter: 1.5mlx2x4</p>	<p>Adapter: 5mlx4</p>	<p>Adapter: 10mlx4</p>
<p>Adapter: 15mlx4</p>		<p>Adapter: 20mlx4</p>	<p>Adapter: 5ml/7mlx4 Vacuum Tube</p>	



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S.





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>

 **ONE
TREE
PLANTED**
<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument used for observing objects too small to be seen by the naked eye. This is achieved through the use of lenses or a single reflective surface in order to magnify images formed on a close proximity to an object. This is achieved through the use of lenses or a single reflective surface in order to magnify images formed on a close proximity to an object.

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device used to store food, pharmaceuticals, chemicals, biological samples, and other items that require low temperatures. It is one of the most important pieces of equipment in a laboratory. It is used to maintain a controlled environment (refrigerated space), without fluids and substances.

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

