



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

# LABORATORY SECTOR

---

# CENTRIFUGES

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## KALSTEIN® Centrifuges

### Tabletop Low Speed Centrifuge

#### Model YR0141



#### Product Description:

It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.

#### Technical Performance:

- ✓ It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Up to 9 user programs storage meet more demands.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#### Technical Specifications:

<b>Model</b>	YR0141	YR0140-1
<b>Max speed</b>	6000r/min	5000r/min
<b>Max RCF</b>	5120xg	4100xg
<b>Max capacity</b>	16x15ml/24x10ml	
<b>Speed accuracy</b>	± 30r/min	
<b>Time setting</b>	1min~99min59s	
<b>Noise</b>	<65db(A)	
<b>Power supply</b>	AC110V/AC220V, 50Hz-60Hz	
<b>Consumption</b>	250w	
<b>Dimension(L×W×H)</b>	410mm×490mm×310mm	
<b>Packing size(L×W×H)</b>	500mm×580mm×400mm	
<b>Net weight (No rotor)</b>	25kg	



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.





All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467 158849 / +33 690760710 / +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S






**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR0141**



 <p>No.2</p>	<p>Capacity: 16x15ml                  Max Speed: 4200r/min                  Max RCF: 2760xg                  Rotor Angle: ≤90°</p>	4 x4 x15ml	Adaptor: 1.5ml, 5ml/7ml、 10ml
		 <p>Adaptor: 1.5ml 5ml/7ml 10ml</p> <p>4x4x15ml</p>	



#Letsgivingore   
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**



**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR0141**



 <p>No.2</p>	<p>Capacity: 4x50ml          Max Speed:          4200r/min          Max RCF: 2760xg          Rotor Angle: ≤60°</p>	<p>4 x50ml</p>	<p>Adaptor: 4 x2 x1.5ml          4 x5ml/7ml          4 x10ml          4 x15ml          4 x20ml</p>
		 <p>50mlx4</p>	<p>Adaptor: 4x2x1.5ml          4x5ml/7ml          4x10ml          4x15ml          4x20ml</p>



#Letsgivemore ♡  
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>

<https://www.kalstein.fr/>

