



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

# LABORATORY SECTOR

---

# CENTRIFUGES

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





# KALSTEIN® Centrifuges

## Tabletop Low Speed Centrifuge

### Model YR131



#### Product Description:

It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.

#### Technical Performance:

- ✓ It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Up to 9 user programs storage meet more demands.
- ✓ It has self-locking, over-speed detection, overtemperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#### Technical Specifications :

<b>Model</b>	YR131
<b>Max speed</b>	5000r/min
<b>Max RCF</b>	5030xg
<b>Max capacity</b>	4x500ml
<b>Speed accuracy</b>	± 10r/min
<b>Time setting</b>	1min~99min59s
<b>Noise</b>	<65db(A)
<b>Power supply</b>	AC110V/AC220V, 50Hz-60Hz
<b>Consumption</b>	750w
<b>Dimension(L×W×H)</b>	460mm×540mm×340mm
<b>Packing size(L×W×H)</b>	550mm×630mm×430mm
<b>Net weight (No rotor)</b>	31kg



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S



## KALSTEIN® Centrifuges Tabletop Low Speed Centrifuge Model YR131



### Rotor for YR131

	<p>Capacity: 4 x 500ml Max Speed: 4000r/min Max RCF: 3040 xg</p>	<p>500 ml Bucket 500ml PE Bottle</p>	<p>Adaptor: 4 x 18 x 1.5ml, 4 x 19 x 5ml/7ml (Blood Vacuum Tube), 4 x 14 x 10ml (Internal Thread), 4 x 9 x 15ml, 4 x 7 x 20ml, 4 x 100ml, 4 x 250ml</p>
		<p>Square Bucket</p>	<p>Adaptor: 4 x 25 x 1.5ml, 4 x 28 x 5ml (Blood Vacuum Tube), 4 x 20 x 10ml (Internal Thread), 4 x 10 x 15ml, 4 x 4 x 50ml, 4 x 250ml</p>
		<p>4 x 14 x 5ml/7ml Aerosol Cover (Option) 4 x 19 x 5ml/7ml Aerosol Cover (Option) 4 x 24 x 7ml Blood Vacuum Tube Microplate Bucket 4 x 2 x 96 Thickness ≤ 15.5mm x 2</p>	



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.









All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467 158849 / +33 680760710 / +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S




**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR131**



 <p>No.2</p>	<p>Capacity:8x50 ml /4 x100ml          Max Speed:5000r/min          Max RCF:5030xg</p>	<p>4 x100ml</p> 	<p>Adaptor:4x4x1.5ml、 4 x4x5ml/7ml、 4 x2 x5ml、 4 x 3 x 10ml、 4x 15ml、 4 x20ml、 4 x50ml、 4 x100ml</p>  <p>Adaptor: 1.5ml, 5ml 5ml/7ml 4x4x10ml</p>  <p>Adaptor: 1.5ml, 5ml 5ml/7ml 10ml 4x4x15ml</p>  <p>Adaptor: 1.5ml, 5ml 5ml/7ml 10ml 8x4x15ml</p>
		 <p>4x2x50ml</p>	<p>4 x2 x50ml</p>



#Letsgivemore   
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.














All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A.S**




**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR131**



 <p>No.3</p>	<p>Capacity:4 x250ml Max Speed:4000 r/min Max RCF:3200 xg</p>	<p>250 ml Bucket 250ml PE Bottle</p>	<p>Adaptor:4x9x1.5ml、4 x7x5ml、4 x9 x5ml/7ml (Blood Vacuum Tube)、4 x 7 x 10ml (Internal Thread)、4 x 5x 15ml、4 x4 x20ml、4 x50ml、4 x100ml、4 x250ml</p> <div style="display: flex; flex-wrap: wrap; justify-content: space-around;"> <div style="text-align: center;">  <p>250 ml Bucket</p> </div> <div style="text-align: center;">  <p>4x9x1.5ml</p> </div> <div style="text-align: center;">  <p>4x7x5ml</p> </div> <div style="text-align: center;">  <p>4x9x5ml/7ml Blood Vacuum Tube</p> </div> <div style="text-align: center;">  <p>4x7x10ml Internal Thread</p> </div> <div style="text-align: center;">  <p>250ml PE Bottle</p> </div> <div style="text-align: center;">  <p>4x5x15ml</p> </div> <div style="text-align: center;">  <p>4x4x20ml</p> </div> <div style="text-align: center;">  <p>4x50ml</p> </div> <div style="text-align: center;">  <p>4x100ml</p> </div> </div>
--	---	--	---



#Letsgivingore   
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.










All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**



**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR131**



 <p>No.3</p>	<p>Max Speed: 4000r/min Max RCF:3200 xg</p>	<p>Adaptor:4 x12 x7ml、4 x6 x15ml、4 x9 x15ml、4 x50ml、4 x2 x50ml、4 x100ml</p> <div style="display: flex; justify-content: space-around; align-items: flex-start;"> <div style="text-align: center;">  <p>4x12x7ml</p> </div> <div style="text-align: center;">  <p>4x6x15ml</p> </div> <div style="text-align: center;">  <p>4x9x15ml</p> </div> </div> <div style="display: flex; justify-content: space-around; align-items: flex-start; margin-top: 20px;"> <div style="text-align: center;">  <p>4x50ml</p> </div> <div style="text-align: center;">  <p>4x2x50ml</p> </div> <div style="text-align: center;">  <p>4x100ml</p> </div> </div>
--	---	--



#Letsgivingore Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S.**




**KALSTEIN® Centrifuges**  
**Tabletop Low Speed Centrifuge**  
**Model YR131**



 <p><b>No.3</b></p>	<p>Capacity:4 x 18 x 7ml Max Speed:4000r/ min Max RCF:3200xg</p>  <p>4x18x7ml</p>	 <p><b>No.4</b></p>	<p>Capacity:2 x 2 x 96 Max Speed:4000r/min Max RCF:2200xg Microplate Bucket</p>
 <p><b>No.5</b></p>	<p>Capacity:12x 15ml Max Speed:5000r/ min Max RCF:3500xg</p>	 <p><b>No.5</b></p>	<p>Capacity:8 x 50ml Max Speed:5000r/min Max RCF:3500xg</p>



#Letsgivemore   
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

