



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges

Tabletop Low Speed Centrifuge

Model YR133



Product Description:

YR133 Tabletop Low Speed lab Centrifuge is widely used in the field Of clinical medicine, biochemistry, genetic engineering, immunology, etc. It is popular with clinic and Reacher lab.

Product Details:

YR133 Tabletop Low Speed lab Centrifuge is widely used in the field Of clinical medicine, biochemistry, genetic engineering, immunology, etc. It is popular with clinic and reached lab.

Technical Performance:

- ✓ It has ABS housing and LED display.
- ✓ The speed and time can be programable.
- ✓ There are lid safety device, it can make the operation safe.
- ✓ The machine can program memory. It can store the specification Of the last operation.
- ✓ Silastic sealed ring comfort to GMP.US FDA certificated.



Technical Specifications:

Model	YR133 Lab centrifuge
Max speed	4000r/min
Max RCF	2220×g
Max capacity	50ml×6
Speed accuracy	± 30r/min
Time setting	1min~60min
Noise	<65dB(A)
Power supply	110V/60Hz 220V/50Hz
Dimension(LxWxH)	320mm×410mm×320mm
Packing size(LxWxH)	420mm×510mm×420mm
Net weight	12kg



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.







All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges
Tabletop Low Speed Centrifuge
Model YR133



CENTRIFUGES

	<p>NO.1 Capacity:7ml×12 Max Speed:4000r/min RCF:2220×g For 5ml/7ml Vacuum Tube</p>		<p>NO.1 Capacity: 7ml×18 Max Speed:4000r/min RCF:2220×g For 5ml/7ml Vacuum Tube</p>
	<p>NO.1 Capacity:7ml×24 Max Speed:4000r/min RCF:2220×g For 5ml/7ml Vacuum Tube</p>		<p>NO.1 Capacity:10ml/15ml× 12 Max Speed:4000r/min RCF:2220×g Adapter:1.5ml、 5ml/7ml Vacuum Tube</p>



#Letsgivemore ♥
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

