



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges

Tabletop Low Speed Centrifuge

Model YR140-1



Product Description:

It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.

Technical Performance:

- ✓ It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Up to 9 user programs storage meet more demands.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



Technical Specifications:

Model	YR140-1
Max speed	5000r/min
Max RCF	4390xg
Max capacity	32x15ml
Speed accuracy	± 30r/min
Time setting	1s~99min59s/1min~99min
Noise	<65db(A)
Power supply	AC110V/AC220V, 50Hz-60Hz
Consumption	450w
Dimension(L×W×H)	460mm×540mm×340mm
Packing size(L×W×H)	550mm×630mm×430mm
Net weight (No rotor)	28kg



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 69760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Tabletop Low Speed Centrifuge

Model YR140-1



 <p>No.1</p>	<p>Capacity: 4 x 50ml / 4 x 100ml Max Speed: 5000r/min Max RCF: 4390 xg</p>	4 x 100ml	Adaptor: 4x4x1.5ml 、 4 x4x5ml/7ml、 4 x2 x5ml、 4 x 3 x 10ml、 4x 15ml、 4 x20ml、 4 x50ml、 4 x100ml	 4x100ml Adaptor: 4x4x1.5ml 4x2x5ml 4x15ml 4x50ml Adaptor: 4x4x5ml/7ml 4x3x10ml 4x20ml
		4 x 50ml	Adaptor: 4x2x1.5ml 、 4 x5ml/7ml、 4x 15ml、 4 x20ml、 4 x50ml	 4x50ml Adaptor: 4x2x1.5ml 4x5ml 4x15ml Adaptor: 4x5ml/7ml 4x10ml 4x20ml
 <p>No.2</p>	<p>Capacity: 8 x 50ml Max Speed: 4000r/min Max RCF: 2810 xg</p>	4 x 2 x 50ml	Adaptor: 8 x 2 x 1.5ml 、 8 x5ml/7ml、 8 x10ml、 8 x15ml、 8 x20ml	 4x2x50ml Adaptor: 8x2x1.5ml 8x5ml 8x15ml Adaptor: 8x5ml/7ml 8x10ml 8x20ml
 <p>No.2</p>	<p>Max Speed: 4000r/min Max RCF: 2810 xg</p>	4 x 4 x 15ml	Adaptor: 4x4x1.5ml 、 4x4x 5ml/7ml、 4x4x5ml、 4 x4x 10ml、 4 x4x 15ml	
		4 x 6 x 15ml	Adaptor: 4x6x1.5ml 、 4x6x 5ml/7ml、 4x6x5ml、 4x6x 10ml、 4 x6x 15ml	
		4 x 8 x 15ml	Adaptor: 4x8x1.5ml 、 4x8x 5ml/7ml、 4x8x5ml、 4x8x 10ml、 4 x8x 15ml	



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S




KALSTEIN® Centrifuges
Tabletop Low Speed Centrifuge
Model YR140-1



 <p>No.3</p>	<p>Capacity: 4 x 18 x 7ml Max Speed:40 00r/min Max RCF:320 0xg</p>	 <p>No.3</p>	<p>Capacity:4 x 12 x 7ml Max Speed:4000r/ min Max RCF:3200xg</p>	 <p>4x12x7ml</p>
 <p>No.4</p>	<p>Capacity: 2 x 2 x 96 Max Speed:40 00r/min Max RCF:220 0xg Microplat e Bucket</p>			



#Letsgivemore 
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

