



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges

Tabletop Refrigerated Centrifuge


Model YR127



Product Description:

Tabletop refrigerated centrifuge is a routine laboratory centrifuge, widely used in research, education and production departments of biology, chemistry, medicine, etc. It is used in laboratories, scientific research, schools and other departments, and is suitable for rapid separation and synthesis of trace samples.

The tabletop refrigerated centrifuge mainly extracts the required components from the liquid mixture. According to the density of each substance, after high-speed rotation, the dense liquid sinks to the bottom layer, and the small density liquid floats on it, thus stratifying the liquid. Refine the pure things I need. It is used in many fields such as medicine and chemical industry.



Tabletop refrigerated centrifuge technical performance:

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy.
- ✓ Touching panel, programmable operation, The operation data can be set according to demands and stored automatically.
- ✓ LED or LCD display, user-friendly interface, simple and convenient operations.
- ✓ The real-time conversation of Speed and RCF is convenient for operation.
- ✓ It has self-locking device, over-speed safety device, over-temperature safety device, imbalance safety device and automatic alarm device.
- ✓ It has 9 files of acceleration and deceleration level for your choice.
- ✓ Silastic sealed ring comfort to GMP. US FDA certificated.



#Let'sgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Tabletop Refrigerated Centrifuge

Model YR127



Technical Specifications :

Model	YR127
Max speed	18000r/min
Max RCF	23755×g
Max capacity	100ml×4
Speed accuracy	± 30r/min
Time setting	1min~99h 59min /1min~99min
Temperature Setting	-20 ~ +40 °C
Temperature Accuracy	±1 °C
Compressor Unit	Non-CFC R134a
Noise	<65db(A)
Power supply	AC110V/60Hz AC220V/50Hz
Dimension(L×W×H)	600mm×470mm×330mm
Packing size(L×W×H)	720mm×610mm×430mm
Net weight	50kg



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Tabletop Refrigerated Centrifuge
Model YR127



 <p>No.1</p>	<p>NO.1 Capacity: 1.5ml/2.2ml×12 Speed:18000r/min RCF:22530×g Adapters:0.2ml、0.5ml</p>	 <p>No.2</p>	<p>NO.2 Capacity:1.5ml/2.2ml×24 Speed:16000r/min RCF:23755×g Adapters:0.2ml、0.5ml</p>
 <p>No.3</p>	<p>NO.3 Capacity: 5ml×10 Speed: 13000r/min RCF:11400×g Adapters:0.5ml、1.5ml</p>	 <p>No.3</p>	<p>NO.3 Capacity: 5ml×12 Speed: 13000r/min RCF:11400×g Adapters:0.5ml、1.5ml</p>
 <p>No.4</p>	<p>NO.4 Capacity: 2×48 well Speed:3000r/min RCF:1400×g</p>	 <p>No.5</p>	<p>NO.5 Capacity:0.5×36ml Speed:13000r/min RCF:16760×g</p>
 <p>No.5</p>	<p>NO.5 Capacity:0.5×48ml Speed:13000r/min RCF:16760×g</p>	 <p>No.5</p>	<p>NO.5 Capacity:10ml×12 Speed:13000r/min RCF:16760×g Adapters:1.5ml、5ml</p>
 <p>No.5</p>	<p>NO.5 Capacity:15ml×8 Speed:13000r/min RCF:16760×g Adapters:1.5ml、5ml</p>	 <p>No.6</p>	<p>NO.6 Capacity: 50ml×6 Speed:12000r/min RCF:13000×g Adapters:15ml、20ml</p>
 <p>No.6</p>	<p>NO.6 Capacity: 100ml×4 Speed:12000r/min RCF:13000×g Adapters:10ml、20ml</p>	 <p>No.7</p>	<p>NO.7 Capacity: 10ml×4 Speed:6000r/min RCF:3200×g</p>



#Letsgivemore ♥
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

