



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR WATER BATHS AND CIRCULATORS



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.



Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.



Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...



#Letsgivemore ♡

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.



Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.



Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.



AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.



KALSTEIN® Water Baths and Circulators Thermostatic Magnetic Stirring Bath

Model YR02311 // YR02312



Product Description:

- ✓ It can be used as water bath or oil bath.

Features:

- ✓ Stainless steel 304 heater.
- ✓ The magnetic stirring system drives the stirrer to rotate synchronistic ally, so that the solution in the bath can be heated and stirred evenly.
- ✓ DC brushless motor ensures stable performance and continuous working.
- ✓ High temperature magnet can continuously work at 300°C without losing magnetism.
- ✓ PID temperature controller ensures accurate temperature control.
- ✓ Key setting and digital display make it easy to operate.
- ✓ YR02311 // YR02312 are equipped with two sensors make bath temperature and container temperature display alternately

Model	YR02311	YR02312
Temperature Range (°C)	Room Temperature+5 ~ 200	
Temperature Stability (°C)	±1	
Temperature Display	Keypad Input, Digital Display	
Stirring Speed Setting	Knob Setting	
Stirring Speed (rpm)	0~2000	
Bath Dimensions (mm)	Φ 220*110	Φ 254*130
Bath Capacity (L)	4	6.5
Heating Power (W)	500	1050
Max. Flask can be placed (ml)	3000	5000
Power Supply (V/Hz)	1~, 11 0V,60Hz or 220-240V~,50/60Hz	
Dimensions (mm)	260W*280D*230H	260W*280D*260H
Net Weight (kg)	5	6



#Let'sgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Water Baths and Circulators Thermostatic Magnetic Stirring Bath

Model YR02311 // YR02312



Product Description:

- ✓ It can be used as water bath or oil bath.

Features:

- ✓ Stainless steel 304 heater.
- ✓ The magnetic stirring system drives the stirrer to rotate synchronistic ally, so that the solution in the bath can be heated and stirred evenly.
- ✓ DC brushless motor ensures stable performance and continuous working.
- ✓ High temperature magnet can continuously work at 300°C without losing magnetism.
- ✓ PID temperature controller ensures accurate temperature control.
- ✓ Key setting and digital display make it easy to operate.
- ✓ YR02311 // YR02312 are equipped with two sensors make bath temperature and container temperature display alternately



#Letsgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF KALSTEIN
EQUIPMENT**

TO WHICH YOU CONTRIBUTE:



**We present the
prototype of the
Kalstein phone**

Stay ahead of
trends with
our advanced
mobile
technology.



www.kalstein.net

www.kalstein.net/en/



PARIS - FRANCE

HEADQUARTERS

5 Rue Castiglione
75001 Paris - France

Tlf: +33 (0) 1 70 39 26 50

Cel: +33 (0) 6 80 76 07 10

Cel: +33 (0) 6 63 81 00 23

sales@kalstein.eu

<https://www.kalstein.eu/>

<https://www.kalstein.net/>

Follow us on our Social Networks!



[@kalsteinlatinoamerica](https://www.instagram.com/kalsteinlatinoamerica)



[@kalsteinlab](https://www.twitter.com/kalsteinlab)



[@kalsteinlabfr](https://www.instagram.com/kalsteinlabfr)

Here you will find our Tutorials!



[@kalsteinenglish4897](https://www.youtube.com/channel/UCkALsteinenglish4897)

