



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR
PATHOLOGICAL
ANATOMY



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.



Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.



Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.



Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.



Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.



AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.



KALSTEIN® Pathological Anatomy Tissue Embedding & Cooling System

Model YR441



Product Description:

- ✓ Fully programmable computer controls allow automatic system start and stop anytime (weekly)
- ✓ Temperature is controlled by microprocessors made in the USA and they are displayed using color-changing LEDs to enable clear visibility of working status
- ✓ Five heated areas, including Paraffin Chamber, Paraffin Dispenser, left and right Thermal Storage Compartments, and heating plate (working area), are individually controlled and work independently without interference from each other
- ✓ Flexible heating mechanism overcomes the shortcomings of traditional technology that can result in excessive temperature differences. System provides fast heating and precise temperature control. In addition, the dual-protection from overheating is safe, reliable and energy-saving



#Letsgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



**KALSTEIN® Pathological Anatomy
Tissue Embedding & Cooling System**

Model YR441



Technical Description:

Model	YR441
Paraffin Chamber Capacity	3 liters
Temperature Range of Forceps Wells	55 – 70°C
Temperature Range of Paraffin-melting Chamber	55 – 70°C
Temperature Range of Thermal Storage Compartments	55 – 70°C
Temperature Range of Heated Working Areas	55 – 70°C
Temperature Range of Paraffin dispenser	55 – 70°C
Temperature Control Precision	±1%
Paraffin Flow Control	Paraffin flow control via finger touch plate and optional foot pedal
Working Temperature of Cryo-Module	≤ -20°C
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	650W / 300W



#Let'sgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF KALSTEIN
EQUIPMENT**

TO WHICH YOU CONTRIBUTE:




**We present the
prototype of the
Kalstein phone**

Stay ahead of
trends with
our advanced
mobile
technology.



<https://onetreepanted.org>



www.kalstein.net

www.kalstein.net/en/



PARIS - FRANCE

HEADQUARTERS

5 Rue Castiglione
75001 Paris - France

Tlf: +33 (0) 1 70 39 26 50

Cel: +33 (0) 6 80 76 07 10

Cel: +33 (0) 6 63 81 00 23

sales@kalstein.eu

<https://www.kalstein.eu/>

<https://www.kalstein.net/>

Follow us on our Social Networks!



[@kalsteinlatinoamerica](https://www.instagram.com/kalsteinlatinoamerica)



[@kalsteinlab](https://www.twitter.com/kalsteinlab)



[@kalsteinlabfr](https://www.instagram.com/kalsteinlabfr)

Here you will find our Tutorials!



[@kalsteinenglish4897](https://www.youtube.com/channel/UCkALsteinenglish4897)

