



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

VETERINARY SECTOR

VETERINARY MONITOR

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Veterinary Monitor Veterinary Blood Pressure Monitor Model YR06110



Product Description:

In view of the pain points such as weak pulse of small animals and inaccurate measurement, animal trembling and agitation leading to measurement failure, shaving is required to accurately measure animal blood pressure, and single point measurement cannot form trend records, Kalstein independently developed the animal electronic blood pressure monitor YR06110, which can easily and quickly test the blood pressure of different sizes of animals, without anesthetic or shaving, to avoid pet shock. Animals quickly enter the test state, one-click operation, intelligent silent pressure, to provide veterinarians with convenient and efficient blood pressure detection equipment.



Features:

- ✓ Accurate and reliable: Unique motion tolerance technology, step-down measurement, anti-jitter interference
- ✓ Large and small animals: Automatically distinguish large and small animals according to the pet's weight
- ✓ Multiple modes: single, continuous, 2 minutes/time multiple measurement modes
- ✓ Comfortable and durable: The soft TPU cuff is more comfortable and sensitive than the traditional rubber cuff
- ✓ Mute measurement: Intelligent mute pressure, broadcast can be muffled, to avoid pets frightened
- ✓ Multiple languages: Chinese, English, Russian, etc.
- ✓ Comprehensive monitoring: pulse, systolic blood pressure, diastolic blood pressure, average blood pressure, trend chart, comprehensive understanding of pet indicators
- ✓ Long standby: Large capacity battery, long standby
- ✓ Easy to carry: built-in lithium battery, no external power supply, easy to move during the test



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Veterinary Montitor Veterinary Blood Pressure Monitor Model YR06110



Technical Description:

Model	YR06110
Display	4.3-inch TFT screen
Measurement method	oscilloscope method
Measurement mode	single measurement, continuous monitoring, 2 minute interval measurement
Blood pressure measurement range	0~280mmHg
Pulse measurement range	0~300 times/min
Blood pressure measurement accuracy	±3 mmHg
Pulse measurement accuracy	±5%
Animal cuff	5
Weight	1.4kg
Size	18cm*15cm*17cm
Power Supply	DC 9.0V power adapter



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

