



Kalstein®

A different accompaniment, at your service

All rights reserved © KALSTEIN France S. A. S.,

VETERINARY SECTOR

VETERINARY PUMPS



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Veterinary Pumps

Veterinary Infusion Pump

Model YR06183



Product Description:

- ✓ 3.5-inch touch screen easy to learn interface, short learning curve.
- ✓ YR06183 satisfies both drug and blood transfusion to meet more clinical purposes.
- ✓ The large-capacity battery of 12 hours at 25ml/h ensures safe transfer.
- ✓ Optional warming strip brings one more option for heating and better heating effect.
- ✓ Optional hook: support hanging on the animal cage for better space management.
- ✓ Highly-advanced ingress protection IP44 effectively blocks water and dust.



Model	YR06183
Infuse mode	Rate mode , Time mode, Weight mode, Intermittent mode, Drug library,Drip mode
Applicable infusion set	Various brands of standard infusion sets
Accuracy	±5%
Infusion rate	(0.10-2000.00)ml/h, increment at 0.01ml/h
VTBI	(0.10~9999.99)ml
Volume infused	(0.00~99999.99)ml
KVO rate	(0.10~5.00ml/h in increments of 0.01ml/h
BOLUS rate	(0.10-2000.00)ml/h
Purge rate	(0.10-2000.00)ml/h in increments of 0.01ml/h
Bubble detection	Ultrasound sensor detection with 6 levels: 25ul-800ul
Occlusion pressure	13 levels
Drug library	Applicable for each infusion mode,and 5000 drugs can be added
History record	50000
Alarm	Air bubble , Occlusion , Door open, Battery exhaust, Low battery , No tube, Finished,Almost done, No operate,AC fail,Empty , Drop abnormal, Standby end , Tube error,Battery error,Overheating, Heating bar not installed , System error. etc.



#Letsgivingmore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most essential equipment that allows to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

