




Kalstein®

A different accompaniment, at your service

All rights reserved © KALSTEIN France S. A. S.,

VETERINARY SECTOR

VETERINARY ELECTROSURGICAL UNIT



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Veterinary Electrosurgical Unit
Veterinary Portable Surgery Diode Laser Systems
Model YR06133 // YR06134



Product Description:

The new portable surgery diode laser system is the best value on the market today. Small size body case with 15W power this mini surgery laser system has 810nm and 980nm dual-wavelengths, 810nm/980nm with pilot beam 300mW for choice, it can be used in the continuous as well as pulsed mode with a contact or non-contact handpiece.

Features:

- ✓ Compact design makes it small in size and light in weight, easy to transport;
- ✓ Highly integrated and modular design;
- ✓ New design inside to make maintenance easy and low cost;
- ✓ Advanced software with presetting protocols to be manually customized by operator;
- ✓ Built-in power meter can test the power in the software by user;
- ✓ Output power calibration at the end of fiber to make the maintenance easy and low cost;
- ✓ Adopt air cooling system guarantee stable performance and long life time.



Technical Description:

Model	YR06133	YR06134
Wavelength	810nm/980nm	810nm+980nm
Maximum Power	15W	15W+15W
Operation Mode	CW, Single or Repeat Pulse	
Pulse Duration	25uS-10S	
Repetition Rate	0.05Hz-20KHz	
Pilot Beam	Red Diode Laser of 635nm, Power	
Control Mode	Touch Screen, 8" True Color	
Dimension	315(H)*215(W)*245(L)mm	
Life time	Over 10000 HRS	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

