



Un accompagnement différent à votre service



Kalstein®

Tous droits réservés © KALSTEIN France S. A. S

LABORATOIRE SECTEUR CENTRIFUGES

www.kalstein.net/fr/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Avantages et soutien

Chez Kalstein France, nous avons à cœur de satisfaire pleinement nos clients, c'est pourquoi nous fournissons des services à valeur ajoutée du plus haut niveau basés sur notre expertise..



Inductions et formations en ligne

Où que vous soyez dans le monde, recevez votre initiation ou votre formation de la part de notre équipe d'ingénieurs dévoués.



Réponse rapide

Notre équipe est toujours disponible pour répondre à toutes vos questions, afin de vous aider dans toutes les situations.



#Letsgivemore ♥

Grâce à votre achat, un don sera fait à une fondation à but non lucratif qui lutte contre le cancer du sein et aide les communautés les plus vulnérables.



Support technique

Bénéficiez de conseils personnalisés pour une maintenance préventive et corrective de vos équipements, grâce à nos manuels et articles Kalstein, aux catalogues spéciaux et aux tutoriels vidéo.



Logistique d'expédition

Nous prenons en charge toute la logistique nécessaire à l'expédition de vos produits, que ce soit par voie maritime, terrestre ou aérienne.



Kalstein Worldwide

Avec plus de 25 ans de croissance avec nos clients, le contenu moderne et multiformat de Kalstein est aujourd'hui présent dans plus de 10 pays.



KALSTEIN® Centrifuges**Centrifugeuse d'étable à grande vitesse****Modelé YR137-3****Description du produit:**

Nous sommes une entreprise professionnelle engagée dans la fabrication et la vente de la centrifugeuse de table à haute capacité YR137-3. Sa conception unique et sa qualité stable lui valent les éloges des clients dans les domaines des collèges et universités, des unités de recherche, des produits pharmaceutiques, des traitements médicaux et de l'ingénierie pétrochimique.

Nous sommes une entreprise professionnelle engagée dans la fabrication et la vente de la centrifugeuse de table de grande capacité YR137-3. Sa conception unique et sa qualité stable lui valent les éloges des clients dans les domaines des collèges et universités, des unités de recherche, des produits pharmaceutiques, des traitements médicaux et de l'ingénierie pétrochimique.

**Performances techniques :**

- ✓ Commande par micro-ordinateur, moteur sans balai à courant continu, fonctionnement stable, faible bruit et précision de la vitesse.
- ✓ panneau tactile, opération programmable, les données d'opération peuvent être réglées en fonction de la demande et stockées automatiquement.
- ✓ Affichage LED ou LCD. Interface conviviale, opérations simples et pratiques.
- ✓ La conversation en temps réel de la vitesse et de la FCR est pratique pour le fonctionnement.
- ✓ L'appareil est doté d'un dispositif de verrouillage automatique, d'un dispositif de sécurité du couvercle, d'un dispositif de sécurité en cas de survitesse et d'un dispositif d'alarme automatique, etc.
- ✓ Il y a 10 fichiers de levier d'accélération et de décélération pour vous permettre de choisir.
- ✓ L'anneau scellé en silicone est conforme aux BPF et certifié par la FDA des États-Unis..

**#Let'sgive more** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Centrifugeuse d'étable à grande vitesse
Modelé YR137-3



Spécifications techniques:

Modèle	YR137-3
Vitesse maximale	5500r/min
FCR max.	5310×g
Capacité maximale	500 ml×4
Précision de la vitesse	± 30r/min
Réglage de l'heure	1min~99min/1min~99h59min
Bruit	<65 dB(A)
Alimentation	110 V/60 Hz 220 V/50 Hz
Dimension (LxlxH)	460mm×540mm×340mm
Taille d'emballage (L x l x H)	550mm×630mm×430mm
Poids net	31kg



#Letsgivemore ♥
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S
















KALSTEIN® Centrifuges


Centrifugeuse d'étable à grande vitesse

Modelé YR137-3



<p>No.4</p>  <p>Rotor capacity: 100mlx4 Max Speed: 5000r/min Max RCF: 4390xg</p>	 <p>100mlx4 Bucket</p>	<p>Adapters: 1.5mlx4x4</p>	<p>Adapters: 5ml/7mlx4x4 Blood Vacuum Tube</p>		
		<p>Adapters: 5mlx2x4</p>	<p>Adapters: 10mlx3x4 G RIA Tube</p>		
		<p>Adapters: 15mlx4</p>	<p>Adapters: 20mlx4</p>		
		<p>Adapters: 50mlx4</p>			
<p>No.5</p>  <p>Rotor capacity: 4200r/min Max RCF: 3100xg</p>	 <p>50mlx2x4 Bucket</p>	<p>Adapters: 1.5mlx2 5ml/7ml Blood Vacuum Tube 5ml, 10ml 15ml, 20ml</p>	 <p>15/10mlx6x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	
		 <p>15/10mlx4x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	 <p>15/10mlx8x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>
<p>No.6</p>  <p>Rotor capacity: 4000r/min Max Speed: 3310xg Max RCF: 500mlx4</p>	 <p>500ml Conical Tube 500mlx4</p>	 <p>500mlx4 V-bottom</p>	 <p>250mlx4 V-bottom</p>		
		 <p>50mlx4x4 V-bottom</p>	 <p>15mlx7x4 V-bottom</p>		



#Letsgivemore 
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges


Centrifugeuse d'étable à grande vitesse

Modelé YR137-3



<p>No.2</p>  <p>Rotor capacity: 500mlx4 Max Speed: 4000r/min Max RCF: 3220xg</p>	 500ml Bucket	 1.5mlx20x4	 5ml/7mlx20x4 Blood Vacuum Tube	 10mlx15x4 GRIA Tube
	 500ml Bottle	 15mlx8x4	 20mlx9x4	 50mlx4x4
	 Bucket	 100mlx2x4	 250mlx4	
	 4x2x96 Microplate	 1.5mlx25x4	 5mlx28x4 Blood Vacuum Tube	 10mlx20x6
	 4x24x7ml Blood Vacuum Tube	 15mlx10x4	 50mlx4x4	 250mlx4



#Letsgivemore 
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Centrifugeuse d'étable à grande vitesse

Modelé YR137-3



<p>No.4</p>  <p>Rotor capacity: 100mlx4 Max Speed: 5000r/min Max RCF: 4390xg</p>	 <p>100mlx4 Bucket</p>	<p>Adapters: 1.5mlx4x4</p>	<p>Adapters: 5ml/7mlx4x4 Blood Vacuum Tube</p>		
		<p>Adapters: 5mlx2x4</p>	<p>Adapters: 10mlx3x4 GRIA Tube</p>		
		<p>Adapters: 15mlx4</p>	<p>Adapters: 20mlx4</p>		
		<p>Adapters: 50mlx4</p>			
<p>No.5</p>  <p>Rotor capacity: 4200r/min Max RCF: 3100xg</p>	 <p>50mlx2x4 Bucket</p>	<p>Adapters: 1.5mlx2 5ml/7ml Blood Vacuum Tube 5ml, 10ml 15ml, 20ml</p>	 <p>15/10mlx6x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	
		 <p>15/10mlx4x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	 <p>15/10mlx8x4 Bucket</p>	<p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>
<p>No.6</p>  <p>Rotor capacity: 4000r/min Max Speed: 3310xg Max RCF: 500mlx4</p>	 <p>500ml Conical Tube 500mlx4</p>	 <p>500mlx4 V-bottom</p>	 <p>250mlx4 V-bottom</p>		
		 <p>50mlx4x4 V-bottom</p>	 <p>15mlx7x4 V-bottom</p>		



#Letsgivemore ❤️
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





Kalstein

Un accompagnement différent, à votre service

#LETSGIVEMORE 

**AVEC
L'ACQUISITION
D'UN
ÉQUIPEMENT
KALSTEIN**

VOUS APPORTEZ UNE
CONTRIBUTION :



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein

(Réseaux sociaux et sites web officiels)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a common piece of laboratory equipment used for research and education.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device that is one of the most important pieces of laboratory equipment. It is used to store samples in a controlled environment (refrigerated space), without freezing and without drying.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Vidéos de formation Kalstein

(Chaîne Youtube et sites web officiels)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

